



Camera Suitable for Trinocular microscopes



Features

➤ **Camera:** Cilika View 12 MP, Sony IMX676

- USB3.0,
- Pixel size: 2μm
- Sensor size: 1/1.6"
- Frame rate: 28fps@3536 × 3536,
- C-mount will be provided
- Cilika View Software: Programmed interval captures (time lapse), video capture by time setting, easy measurement calibration, Measurement in microns, inches, and mm, length measurement, ellipse, rectangular, irregular shape measurement, Perimeter, radius, circumference, angle measurement, zoom function, online file sending, and receiving data export, report generation and printout, extended depth of field (EDF).

Manufactured by:

MedPrime Technologies Pvt. Ltd.

Address: 2nd Floor, Plot no. A-246, Road no. 30-B, Wagle Estate, Thane, Maharashtra 400604

Website: www.medprimetech.com

Phone: +91 8976924214